



FREE  
DOWNLOAD

Google play



twion

**twion**<sup>®</sup>  
Alber Smartphone App



Simply mobile.

twion Mobility App: Help for underway and additional features on the twion

## THE CLEVER HELP FOR UNDERWAY

It provides you with valuable tips and tricks for daily life with your Alber product. In case of a breakdown, identify the cause of the fault using the App and alleviate it yourself in many cases. If that isn't possible for some reason, contact the Alber Service Centre directly from the App or your local specialist dealer. In this way, you always get help fast, directly on your smartphone – no matter where you are at that moment.

### Functions

- Quickstart guide, brief operating instructions
- List of specialist dealers worldwide
- Troubleshooting with solution management
- Contact specialist dealer or Alber Service Centre directly
- Practical tips on all aspects of the twion



## MOBILITY PLUS PACKAGE

With the Mobility Plus package you can activate and enjoy even more features in the twion Mobility App and communicate with the twion wheels via smartphone!

Prerequisite is a Smartphone with Android operating system and ability to connect to the internet.

The Mobility Plus package contains an outer ring with an exclusive design, a booklet with operating instructions, and a licence code to activate the other features.



FREE  
DOWNLOAD



## AN OVERVIEW OF THE FUNCTIONS



### Driving information

- Readout of running parameters including battery status,  $\emptyset$  speed, maximum speed, current speed, altimeter, odometer
- Simple switchover between ECO (less motor assistance for greater range) and SPORT (more dynamic ride characteristics with longer run-on) modes
- Increase assist up to 10 km/h



### Silent mode

- Mute mode: by simply pressing a button you can turn off the charge indicator and the acoustic signals from the wheels
- An individual profile can be set up



### Tour mode

- Recording and display of journeys using a GPS signal
- Detailed display of TOUR data such as TOUR distance covered, Tour time, altimeter and number of calories burned
- Tour data is stored in chronological order



### BlueDrive controller

- When the wheelchair is unoccupied it can be controlled remotely using a smartphone
- You can choose whether to control it with the built-in inclination sensor (by gently tilting the smartphone) or with a virtual joystick



EXPAND YOUR TWION APP WITH MOBILITY PLUS PACKAGE AND ENJOY MORE FEATURES



Alber GmbH  
Vor dem Weißen Stein 21  
72461 Albstadt  
Germany  
Phone +49 7432 2006-0  
Fax +49 7432 2006-298  
[www.alber.de](http://www.alber.de)

