

- Standard
 Mandatory fields:¹
Information required for processing
 Option field

- ORDER
 QUOTE

Invoice address

cust. no.

PO no.

Company

Street

City State ZIP

Order marked for

Shipping address (if different from invoice address)

Company

Street

City State ZIP

Type of order

Mounting by Dealer² (Please fill out the „Wheelchair details“ section completely!)
 Mounting at Alber³ **\$65.00**

<input type="checkbox"/> e-motion M25 – wheel size 22"	1592520	\$7,733.00
<input type="checkbox"/> e-motion M25 – wheel size 24"	1592521	\$7,733.00
<input type="checkbox"/> e-motion M25 – wheel size 25" <small>Flat Free Tires Not Available for 25" Wheels</small>	1592522	\$7,945.00
<input type="checkbox"/> ECS remote Control	1592486	\$1,193.00

Wheelchair details

Wheelchair manufacturer (e.g. Invacare)

Model (e.g. Pro Spin)

Year of manufacture

Version (e.g. Hemi)


NO DRUM BRAKES – incompatible with Alber Wheels
NO 5TH WHEEL ANTI TIPPERS

Important hardware


Bracket for mounting to the wheelchair 1467178 **\$0.00**
 Anti-tippers with jack-up function 1593206 (pair) **\$633.00**

Existing wheel position


1 = back
 2 = center
 3 = front



Seat width: (seat tube measured externally) in



Seat height back cm (measured to lower edge of seat sling) in



Wheel camber (if possible)

0°
 3°

Comments (mounting, delivery etc.)

- Standard **Mandatory fields:**¹
Information required for processing
 Option field

Tires

- Flat Free Tire 1 3/8" black 1592533 \$0.00
 Schwalbe Right Run (Standard) black 1592530 \$0.00
(1" for wheel size 22" and 24", 1.1" for wheel size 25")

Push rims

- Stainless steel (standard) 1592526 \$0.00
 Coated push rim Tetra 1592527 \$409.00
(not available for 25")
(+250 g per wheel)
 Carbolife CURVE L, powder coated, black 1592528 \$875.00
(not available for 22")
 Carbolife QUADRO Tetra-Gripp coating 1592529 \$1,034.00
(24" wheel only)

Push rim covers (removable)

- push rim cover 22" 1592425 (pair) \$299.00
 push rim cover 24" 1592421 (pair) \$299.00
 push rim cover 25" 1592426 (pair) \$299.00

Spoke cover

- Spoke cover 22", Design „Comic“ 1592266 (pair) \$347.00
 Spoke cover 24", Design „Sphere“ 1592268 (pair) \$347.00
 Spoke cover 25", Design „Sphere“ 1592269 (pair) \$347.00

Accessories

- Transport Axle 1591362 \$120.00
 Loop Axle 1591500 \$251.00
 Bracket for ECS remote control 1592743 \$265.00
for fitting to the wheelchair incl. adapter for smart phone

Installation of bracket remote (option)

- Right
 Left (instead of right)



Date / /

Signature of ordering party

¹ In order to smoothly and quickly process your order we require ALL information in the fields highlighted in blue.
² Installation kits can only be shipped to trained specialist dealers (installation training).
³ Please send the wheelchair with a copy of this order as soon as possible and freight paid to Alber.

