

PORTABLE POWER JUST GOT BETTER

e-fix Eco

Joystick controlled power
add-on for manual wheelchairs

The e-fix eco offers many of the same features and benefits you can expect from the e-fix product line including:

- Enhanced Mobility and Independence
- Convenience
- Easy Travel
- Attendant Control



HERE'S HOW WE MADE IT AN EVEN BETTER DESIGN

Super quiet motor on the 22" & 24" wheels



New discreet black
battery color



New battery bag style
with straps on pack
helps ensure stability
and strength



Modern dark gray
wheel hub



New flat and modern
design on control unit

PERFORMANCE



- Free-wheeling thanks to the clutch
- Individual care through versatile programming functions built into the control unit
- Control unit is adaptable to be bracketed as the attendant control
- New joystick options provide end-users with more versatile functionality
- Energy recovery braking technology - Energy in the motor transfers back into the battery resulting in long battery life and provides the end-user recovery while pushing

DIFFERENCES BETWEEN THE E-FIX AND THE E-FIX ECO

e-fix

- Comes in two versions. The e-fix (e35) and the e-fix HD (e36) which features a higher max user weight (from 265lbs. to 352lbs.)

e-fix eco

- Comes in one version which features a max user weight of (265lbs.)
- Equipped with a Lead-Acid battery weighing 19.4 lbs. As opposed to the e-fix (e35 and e36) which features a Lithium-Ion battery weighing 4.4 lbs.
- Battery pack must be removed from the wheelchair to fold for transportation as opposed to the e-fix (e35 and e36) which can be folded while the battery remains attached to the wheelchair.



TECHNICAL DATA

	e-fix (e35)	e-fix (e36)	e-fix Eco
Battery	Lithium-ion	Lithium-ion	Lead-acid
Range*	10 miles	10 miles	10 miles
Max User Weight	265 lbs.	352 lbs.	265 lbs.
Max Total Weight	375 lbs.	463 lbs.	375 lbs.
Attendant Control	Yes	Yes	Yes
Plug	Magnetic	Magnetic	Screw
Display	Color	Color	Color
Bluetooth	Yes (EU Market)	Yes (EU Market)	No
Drivewheel Sizes	22" & 24"	24"	14", 22" & 24"
Opti-box	Yes	Yes	No

Component Weights

Drive wheels	17.4 lbs.	17.2 lbs.	17.4 lbs.
Battery pack	4.6 lbs.	4.4 lbs.	19.4 lbs.
Control unit	1.3 lbs.	1.3 lbs.	1.3 lbs.
Interface	1.8 lbs.	1.8 lbs.	1.1 lbs.
Color	White Design	White Design	Gray Design

*Depending on the type of wheelchair, weight of the occupant, type of terrain.
With optional, high capacity battery pack, up to 12.4 miles range

Display - Intuitive and user-friendly color display. Can be individually programmed to meet your driving characteristics.



FOR MORE ABOUT ALBER PRODUCTS:

alber-usa.com/us/ • smoov.com/us-en/

FOR MORE ABOUT INVACARE'S MOBILITY PORTFOLIO:

rehab.invacare.com

©2022 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols™ and®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 21-251 211211



Yes, you can.

